

Washing MinION Flow Cells

Version: WKE_1012_v1_revK_08Apr2016
Last update: 03/10/2017

Flow Cell Number:

DNA Samples:

Before start checklist

☐ Wash Kit (EXP-WSH002)

☐ Pipette tips P1000

☐ Pipettes and pipette tips

☐ Pipettes P1000

INSTRUCTIONS

NOTES/OBSERVATIONS

Starting the wash process

☐ Click the Stop button in the MinKNOW GUI.

☐ Open the priming port cover of the Flow Cell to check that buffer is continuous.

☐ Ensure that the SpotON sample port is closed with the activator, and the priming port is open for the washing steps.

☐ Using a Gilson P1000 or similar add 150 µl of Solution A through the priming port of the Flow Cell.

☐ Wait 10 minutes.

Follow one of the two options described in the next steps of the protocol.

If planning to add the next library immediately

☐ Add 150 µl of Solution B through the priming port of the Flow Cell.

The Flow Cell is now ready to be used as described in the steps from "Loading a Library" onwards in the protocol being used for the experiment.

IMPORTANT

☐ A Platform QC cannot be run on the Flow Cell prior to loading a subsequent library if the above process is used.

If planning to store the Flow Cell for later use

☐ Slowly add 500 µl of Storage Buffer through the priming port of the Flow Cell.

☐ Close the priming port, and remove all buffer from the waste section of the Flow Cell through either of waste ports.

☐ The Flow Cell can now be stored at 4-8 °C.

☐ When you wish to reuse the Flow Cell, remove the Flow Cell from storage, and allow to warm to ambient temperature (~5 minutes).

Carry out Platform QC as described in the "Preparing for an experiment" section of the protocol being used.